

TECHNICAL DATA SHEET

Instrument Pre-Soak Concentrate - Multi-Enzyme

SKU	ACX-ENZPRESOAK-VAR
Brand	Alle's ClinX

1. Product Overview

This publishing-safe TDS provides product identification and non-operational technical information. For handling, hazard classification, PPE, first aid, spill response, storage controls, and any use instructions, refer to the current product label and Safety Data Sheet (SDS).

Product / chemistry	Protease + lipase + amylase enzyme blend
pH	6.5–7.5 (neutral)
Intended sectors	Education, Healthcare, Labs
Shelf life	18 months from date of manufacture, unopened
Pack size	5 Litres / Bulk on request
Water-hardness note	Use softened or RO water where local supply exceeds 300 ppm CaCO ₃ - hardness ions can inhibit enzyme activity

2. Safety Reference

SAFETY: Use only under the current product label, SDS, and an approved site procedure. Do not mix chemical products unless an approved controlled document explicitly permits it.

PPE reference: Gloves recommended when handling soiled instruments

This document intentionally omits dilution, dosing, mixing, dwell-time, and step-by-step use instructions. Those controls belong in the current label/SDS and trained-site procedures.